
WeightWatchers®  
BY Tefal

FR  
NL  
DE  
SV  
DA  
EN

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In association with Weight Watchers, Tefal has developed these bathroom scales specifically designed to track your weight with various indicators including your BMI\*.

\* **BMI = Body Mass Index.**

The BMI measures the weight status of a person. It is used by the World Health Organization (WHO) to assess risks caused by overweight or anorexia in relation to a reference scale.

BMI is calculated as follows:

weight in kilograms / (height in metres x height in metres).

For example, for a person weighing 62 kg and 1.60 metres tall:

$$62 / (1.6 \times 1.6) = 24.2$$

The **BMI is 24.2.**

Tefal and Weight Watchers bring you bathroom scales that automatically calculate your BMI and interpret it to help you track your weight better.

The drawings in this user's manual are intended for illustration only and may differ from the actual product.

## Description

Memory number indicators

Message prompting you to step off the platform

BMI\* LED

Message: low battery – replace them

EN

Clear memory key



OK key

Scroll key (minus)

Scroll key (plus)

Platform: As soon as you step onto the platform, the scales turn on automatically

Display screen:  
– weight  
– height  
– weight difference compared to the previous time  
– BMI

## Tips on tracking your weight

Weight fluctuations can be high during the day, depending on your food intake, digestion or physical activity. We recommend that you **always measure your weight in the same conditions, preferably first thing in the morning before eating any food or drinking any liquid**. There is no need to weigh yourself several times in a single day. Measure your weight preferably every week rather than every day.

## Precautions of use

- Always place the scales on a flat, hard floor when you take your measurement. Avoid carpets or soft surfaces.
- Once you step onto the scales, try to stand still.
- Do not move the appliance when handling the keys.
- Caution: these bathroom scales are calibrated to track the BMI only for people above 18. They are not calibrated for children.
- Caution: all the data stored in the memory will be lost when you replace the batteries.

EN

## Using your bathroom scales

Automatic switch-on when you step onto the platform.

### 1 – Selecting your data

When using for the first time, you must choose your memory number and your height so that the scales can calculate your BMI and your personalised weight difference – **see sequence I, page 3**.

– Step onto the scales.

– Select your memory by pressing , until the chosen memory number is displayed.

Confirm your memory number by pressing the  key for 2 seconds.

– Select your height (100 to 220 cm)

. if you are less than 170 cm tall, press , until your height is displayed.

. if you are more than 170 cm tall, press , until your height is displayed.

Confirm your height by pressing the  key for 2 seconds.

## 2 – Data supplied by your scales

During subsequent uses, the scales will automatically recognize you thanks to its exclusive automatic recognition system and will display your **current weight**, the **weight difference** as compared to the last measurement and your **BMI** – see **sequence II, page 3**.

If 2 users have similar weights or if your weight difference is very high, the scales will ask you to reconfirm your memory number.

For this, after stepping off the platform, press the  $\ominus$  or  $\oplus$  key until your memory number and height are displayed and then confirm by pressing the  $\boxed{\text{ENTER}}$  key for 2 seconds.

### Body Mass Index

- Your scales calculate your BMI automatically and interpret it:

<18: Underweight

18 to 25: Normal weight

25 to 30: Overweight

>30: Obese

## 3 – Clearing the memory

When you want to reset your memory number, you must delete the data contained in it – see **sequence III, page 3**:

- Step onto the scales, which will recognize you.
- Step off the platform, press and hold down the  $\boxed{\text{CLEAR X}}$  key until the “Clr” message is displayed.
- To record your data on the memory, repeat from **sequence I, page 3**.

When all the memory numbers are in use, you must delete the data stored in a memory number to record your data.

- Step onto the scales, which will not recognize you.
- Step off the platform. Press the  $\ominus$  or  $\oplus$  key until the memory number you want to clear is displayed, and then hold the  $\boxed{\text{CLEAR X}}$  key down until the “Clr” message is displayed.
- To record your data on the memory, repeat from **sequence I, page 3**.

## Special messages

Messages/Situations	Actions/Solutions
	The load is higher than the maximum weight of 160 kg. <b>→ Do not exceed the maximum load.</b>
	You must step off the platform to continue the sequence. <b>→ Step off the platform.</b>
	See page 2
	<b>1) All the memory numbers are in use.</b> <b>→</b> see “Clearing the memory”.  <b>2) You are not recognized.</b> <b>→</b> Select your memory number and your height by pressing $\oplus$ , and press the $\boxed{\text{ENTER}}$ key for 2 seconds.

EN

You press a key and there is no effect on the appliance.

**1) The scales are detecting another key than the one pressed by your finger.**

→ Make sure:

- you are pressing the key with just one finger
- you are not resting on the platform
- you are not holding it pinched
- the scales are placed properly on the floor.

Or

**2) You have not pressed the key long enough.**

→ Keep your finger pressed on the key for 2 seconds.

EN

**Environment protection first!**

① Your appliance contains valuable materials which can be recovered or recycled.

➡ Leave it at a local civic waste collection point.