F	p. 2 + 6	Service consommateurs & commande accessoires tél. 0 810 77 47 74 - www.tefal.fr
NL	p. 2 + 10	+ 31 318 58 24 24 - www.tefal.nl
D	p. 2 + 13	+ 49 212 387 400 - www.tefal.de
I	p. 2 + 16	+ 39 199 20 73 54 - www.tefal.it
GB	p. 2 + 19	Consumer Service & Accessories ordering tel. 44 845 725 85 88 - www.tefal.co.uk
E	p. 2 + 22	+ 34 902 312 400 - www.tefal.es
P	p. 2 + 25	Clube Consumidor + 351 808 284 735 - www.tefal.pt
GR	p. 2 + 28	+ 30 801 11 97 101 - www.tefal.com 210 - 6371251
TR	p. 2 + 31	
DK	p. 2 + 34	+ 45 4466 31 55 - www.tefal.dk
S	p. 2 + 37	Kundservice + 46 8 594 213 30 - www.tefal.se
FIN	p. 2 + 40	Asiakaspalvelu + 358 962 29 420 - www.tefal.fi
N	p. 2 + 43	Kundeservice & bestilling av reservedeler + 47 815 09 567 - www.tefal.no

Tefal

BODYUP


* COMPACIE 360 EURO RCCC - 217283 - 707 - RCS TEFLammy B 301501 920 - Ref. 201500349



www.tefal.com

Tefal
www.tefal.com


F
NL
D
I
GB
E
P
GR
TR
DK
S
FIN
N

1st use - 1^{re} utilisation

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The drawings are designed merely to illustrate the features of the scales and do not strictly reflect reality.

In the consumer's interest, TEFAL SAS reserves the right to modify the characteristics or components of the product at any time.

Description

1. Display of the identification numbers of the 4 memories
2. Battery status indicator
3. Selection : Man/woman
4. Cancel key
5. Scroll and delete keys
6. Display of the different values: height in cm (100 to 220 cm), age (starting at 6 years), weight and body fat in kg
7. Graphic representation of body fat and recommended minimum and maximum limits
8. ENTER key
9. Electrodes for positioning the feet

Measuring principle

Caution: this appliance must not be used by people with pacemakers or other medical electronic devices due to the risk of interference from electric current.

When you step on the scales, a very weak electric current is transmitted to your body by the electrodes.

This current, which travels from one foot to the other, is slowed down by body fat. This body resistance, known as bio-electrical impedance, varies according to the user's sex, height and age. It makes it possible to measure body fat and its minimum and maximum range.

This appliance presents no danger as it runs on very weak current (2 LR03 AAA 1.5 V batteries). However, there is a risk of getting a false measurement in the following cases:

- Person with fever
- Child under the age of 6 (appliance not calibrated below that age)
- Pregnant woman (no danger)
- Person under dialysis, person with swollen legs, person suffering from dysmorphia.

Preventing obesity starting in childhood

These scales provide an efficient way of tracking weight problems right from childhood. Indeed, the growth period is when the characteristics of the human body change the most and it is often difficult to detect a drift. Yet a tendency towards obesity, for example, is treated much more efficiently if it is detected early, thereby avoiding that the child becomes obese as an adult. These scales are designed to keep track of the percentage of body fat from the age of 6.

GB


Functions

- The scales provide **customised records**.

They measure and display:

- **Your weight in kg**
- **Your body fat in kg.**

They display these results on a **graph** according to your minimum and maximum range (calculated according to your sex, height and age). The objective is to make sure that your figures are within those limits.

If your body fat is above the upper limit, you should change your eating habits and do physical exercise.

If your body fat is below the lower limit, it would be wise to increase your body fat as it insulates the body, protects it and acts as its energy reserve.

- Your personal characteristics (sex, height, age) are saved in the scales' memory and you are **automatically recognised** by the scales the next time you weigh yourself.

When you change the batteries, all saved data is deleted.

- The scales **can identify up to 4 different users**, including **children under 6**.

Operating advice

These scales must be used in accordance with the **following conditions**:

- **The measurement MUST be taken bare foot** (direct contact with the skin).

The bare feet must be dry and properly positioned on the electrodes.

- **It is important to check that the feet and legs** (thighs, knees, ankles) **do not touch each other**. If need be, slide a sheet of paper between the legs.

- **The body measurement must always be done under the same conditions**, preferably once a week, in the morning (so as not to take account of weight fluctuations throughout the day due to eating, digestion or physical activity) **15 minutes after getting up**, to ensure that the body water has spread to the lower limbs.

- **Do not measure yourself after a physical effort**. Wait 6 or 8 hours to get a significant measurement or measure yourself before the effort.

Modification of personal characteristics (such as children's age or height)

1. When your personal data appears - **fig. 5b**, press for 5 seconds. **enter**
2. The image of the person blinks as well as the word **SELECTION**. If you want to change the sex, press or then press once, otherwise press **enter**
3. The height blinks as well as the word **SELECTION**. If you want to change the height, press or then press once, otherwise press **enter**
4. The age blinks as well as the word **SELECTION**. If you want to change the age, press or then press once, otherwise press **enter**

The results of your previous weigh-in will be displayed automatically then the scales will switch off.

GB

Environment protection first!

- Your appliance contains valuable materials which can be recovered or recycled.
- Leave it at a local civic waste collection point.

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