# Tefal®

## GOAL™

EN


RU

UK


2015900577


[www.tefal.com](http://www.tefal.com)

## LR03 AAA

I

Goal

II



Bip





Bip

Bip

**EN.** The drawings in this user's manual are intended for illustration only and may slightly differ from the product.

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This bathroom scale automatically calculates your BMI\* and brings you a clear read so that you can better track your weight.

### \*BMI = Body Mass Index

The BMI measures the weight status of a person. It is used by the World Health Organization (WHO) to assess risks caused by overweight or anorexia in relation to a reference scale.


#### **BMI is calculated as follows:**

Weight in kilograms / (height in metres x height in metres).

For example, for a person weighing 62 kg and 1.60 metre tall:

**62 / (1.6 x 1.6) = 24.2 The BMI is 24.2.**

## Description

## Tips on tracking your weight

Weight fluctuations can be high during the day, depending on your food intake, digestion or physical activity. We recommend that you always measure your weight in the same conditions, preferably first thing in the morning before eating any food or drinking any liquid. There is no need to weigh yourself several times in a single day. Measure your weight preferably every week rather than every day.

## Precautions of use

- **Always place the scale on a flat & hard floor when you take your measurement. Avoid carpets and soft surfaces.**
- **Once you step onto the scale, try to stand still.**
- **Do not move the appliance when handling the keys.**
- **Caution: this bathroom scale is calibrated to track the BMI only for people above 18. It is not calibrated for children.**
- **Caution: all the data stored in the memories will be lost when replacing batteries.**

## Using your bathroom scale

Automatic switch-on when stepping onto the platform.

### 1. Selecting your data

When you measure your weight for the first time, you must set your memory number, height and goal weight, so that the appliance can calculate your BMI and the progress to make to achieve the weight you have set as your goal – **see sequence I, page 3.**

- Step onto the bathroom scale.
- Select your memory by pressing  until the chosen memory number is displayed. Confirm your memory number by pressing the  key for 2 seconds.
- Select your height (100 to 220 cm):
  - if you are less than 170 cm tall, press  , until your height is displayed.
  - if you are more than 170 cm tall, press  , until your height is displayed.Confirm your height by pressing the  key for 2 seconds.

- Select your goal weight:
    - If you want to achieve a weight lower than your current weight, press  until you reach your goal weight.
    - If you want to achieve a weight higher than your current weight, press  until you reach your goal weight.
- Confirm your goal by pressing the  key for 2 seconds.

The illuminated progress bar is activated.

**Caution:** the goal weight should be reasonable. You cannot record a goal weight corresponding to a BMI of less than 18, considering your height.

## 2. Data supplied by the appliance

The appliance is equipped with an automatic user recognition system.

During subsequent uses, the appliance will automatically recognize you thanks to its exclusive automatic recognition system and will display your current weight, the progress made (as a value) in achieving your goal weight compared to the last measurement and your BMI – **see sequence II, page 3.**

If 2 users have similar weights or if your weight difference is very high, the appliance will ask you to reconfirm your memory number.

For this, step off the platform, press the  or  key until your memory number and height are displayed and then confirm by pressing the  key for 2 seconds.

## Body Mass Index

- The appliance calculates your BMI automatically and interprets it:

<18: Underweight

18 to 25: Normal weight

25 to 30: Overweight

>30: Obesity

## Tracking the goal weight

Every time you measure your weight, depending on the selected goal, the progress made in achieving the stored goal weight is displayed as a value on the LCD screen and on the illuminated progress bar.

For example: ***You have a weight loss goal***

## You have a weight gain goal

**Important: The progress made** **is always calculated from your weight at the time you selected your goal weight.**

In case you do not want to enter this goal weight, the difference in weight from your previous weighing is displayed on every new measurement, and the illuminated progress bar remains inactive.

### Modification or cancellation of the goal weight

#### - Cancellation

When your goal weight or the illuminated progress bar is displayed, press for 10 seconds until you hear 3 beep sounds.

#### - Modification

When your goal weight or the illuminated progress bar is displayed, press for 2 seconds until the goal weight flashes. Refer to

**“Select your goal weight” page 9.**

### 3. Memory reset

When you want to reset memory to change your height details, you must delete the data contained in the memory – **see sequence III, page 3:**

-Step onto the bathroom scale, you are recognized.

-Step off the platform, press and hold down the  key until the "Clr" message is displayed.

-To record your data, repeat **from sequence I, page 3.**

When all the memory numbers are in use, you must delete the data stored in a memory to record your data.

-Step onto the bathroom scale, you are not recognized.

-Step off the platform. Display the memory number you want to clear by pressing the  or  key, and then hold the  key down until the "Clr" message is displayed.

-To record your data on the memory, repeat from **sequence I, page 3.**

## Special Messages

Messages / Situations	Actions / Solutions
	The load exceeds the maximum capacity of 160 kg. » <b>Do not exceed the maximum capacity.</b>
	The appliance needs you to step off the platform to continue the sequence. » <b>Step off the platform.</b>
	See page 2

**1) All the memories are used.**

» See "memory reset".

**2) The bathroom scale does not recognise you.**

» Select your memory and your height by pressing on  then press 2 seconds on  .

You press a button to save your data and the scale does not react.

**1) The appliance detects something other than your finger pressing.**

» Ensure:

- that you are pressing with just one finger,
- that you are not leaning on the platform,
- that you are not holding it,
- that the appliance is placed on the ground.

or

**2) You have pressed the button too briefly.**

» Keep your finger pressing 2 seconds on the button.

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